# Enbio S Enbio PRO

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User manual EN 31

**CE** 2274  
ver. 26.07.2021

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





RYS. 3. REAR VIEW.

## Table of Contents

<b>1. Introduction .....</b>	<b>33</b>	6.2 Periodic inspections.....	48
1.1 Purpose.....	33	<b>7. Data archiving.....</b>	<b>49</b>
1.2 Applicable legal acts.....	33	<b>8. ENBIODATAVIEWER software .....</b>	<b>49</b>
1.3 Purpose of the device.....	33	8.1 Software installation.....	50
1.4 Symbols used on the device.....	33	8.2 Program construction and main functionalities.....	52
1.5 Precautions, requirements and recommendations.....	34	<b>9. Warning messages and error codes.....</b>	<b>55</b>
<b>2. Scope of delivery and unpacking of the device .....</b>	<b>35</b>	9.1 Warning messages.....	55
2.1 Unpacking of the device.....	35	9.2 Information Codes.....	55
2.2 Standard equipment.....	35	9.3 Error codes.....	56
<b>3. Device installation.....</b>	<b>35</b>	<b>10. Warranty claim handling .....</b>	<b>57</b>
3.1 Water quality.....	36	<b>11. Warranty terms and conditions .....</b>	<b>58</b>
<b>4. Tool preparation and loading .....</b>	<b>37</b>	<b>12. Technical information.....</b>	<b>58</b>
4.1 Tool pack preparation.....	37	<b>13. EC Declaration of conformity.....</b>	<b>59</b>
4.1.1 Characteristics of a sterilisation pack.....	37		
4.1.2 Rules for arranging tools on a tray .....	37		
4.1.3 Principles of packing tools for sterilisation.....	38		
<b>5. Starting the device.....</b>	<b>38</b>		
5.1 Commissioning.....	38		
5.2 Program selection.....	38		
5.3 Test programs.....	41		
5.4 Information menu.....	44		
5.4.1 Counters.....	45		
5.5 Restarting .....	45		
<b>6. Maintenance and care .....</b>	<b>47</b>		
6.1 Replacement parts.....	48		

The latest version of the manual is available at [www.enbio.com](http://www.enbio.com)

# 1. INTRODUCTION

## 1.1 Purpose

The purpose of this user manual is to provide information about the ENBIO S/ ENBIO PRO steriliser and ensure:

- proper installation and setup,
- optimum use,
- safe and reliable operation,
- regular and correct maintenance and servicing in accordance with requirements.

## 1.2 Applicable legal acts

The ENBIO S/ ENBIO PRO sterilisers meet the following legal requirements:

- The sterilisers have been designed and manufactured in accordance with the EN 13060 standard.
- The sterilisers meet the essential requirements of Directive 93/42/EEC and the Medical Products act and have the status of medical product.
- According to Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE)

## 1.3 Purpose of the device

ENBIO S and ENBIO PRO sterilizers carry out sterilization processes in class B and class S.

ENBIO S and ENBIO PRO devices are small class B steam sterilizers in accordance with EN 13060, classified as a medical device of class IIb, in accordance with Annex IX of Directive 93/42 / EEC concerning medical devices and the Regulation of the Minister of Health of November 5, 2010. on how to classify medical devices. According to the 13060 classification, medical devices can be sterilized in ENBIO S and ENBIO PRO devices: solid loads, small porous objects, small porous loads, full porous loads, simple objects recessed, narrow-clearance items, multiple packaging that can be unwrapped or wrapped (single and multi-layered).

**The 134 FAST process is dedicated to solid, non-porous, simple instruments and dental tools (eg scissors, handles, pliers, chisels, probes, etc.) exclusively unpacked, not textile.**

In addition, the Enbio PRO autoclave has a dedicated PRION program as one of the stages of decontamination of items that have or may have had contact with diseased prion proteins (e.g. Creutzfeld-Jacob disease, BSE, etc.): solid loads, small porous items, small porous ripples, full porous loads, simple hollow items, items with narrow clearance, multiple packages that may be unpacked or wrapped (single and multi-layer).

The ENBIO S and ENBIO PRO devices can be used in primary healthcare practices, in dental practices, in treatment rooms. The sterilizer is adapted to work in the vicinity of other electrical medical devices.

Non-medical application:

The ENBIO S and ENBIO PRO devices can also be used for non-medical applications, such as in beauty and wellness salons, as well as in veterinary practices, tattoo studios, earrings and hairdressing salons.

 **The ENBIO S/ ENBIO PRO devices can't be used to sterilise liquids, biomedical waste or pharmaceutical products.**

The device is intended for professional use by properly trained staff only.

## 1.4 Symbols used on the device

This symbol is located on the front of the device, on the upper part of the drawer front. It is recommended to maintain caution due to high temperature within and around the operating chamber.

This symbol is located on the device's rating plate and indicates the serial number.

This symbol is located on the device's rating plate and indicates compliance with EC guidelines.

This symbol is located on the device's rating plate and indicates the device's date of production.

This symbol is located on the device's rating plate and indicates the device's manufacturer.

This symbol is located in the user manual and indicates reading the information provided in the user manual.

According to Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE) it is not allowed to place, dispose of, store used electrical and electronic equipment along with other waste. Used equipment should be delivered to a local collection point for used electrical equipment, which is registered in the Chief Inspectorate for Environmental Protection and conducts selective waste collection.

### 1.5 Precautions, requirements and recommendations

- The user is responsible for the installation, correct operation and maintenance of the device in accordance with instructions provided in this user manual. If needed, contact the service or the supplier of the device.
- The steriliser is not intended for sterilising liquids, biomedical waste or pharmaceutical products.
- The steriliser must not be used if explosive gases or vapours are present in the air.
- After the cycle is completed, the load is hot. Remove tools or packs from the chamber using appropriate thermal gloves or equipment that prevents burns.
- Do not remove the rating plate or any other elements of labelling from

the device.

- Follow guidelines for preparing tools for sterilisation.
- Pouring water or other liquids on the device may cause a short-circuit.
- Prior to inspection, maintenance or servicing, turn off the device and disconnect it from the power grid.
- Servicing may only be performed by trained service personnel and using original replacement parts.

**After sterilization tools are hot and wet in the 134 FAST program.**

Read this Operating Manual carefully before using this device. Install and operate the device strictly as specified herein. Comply with all safety requirements for the device. This will ensure proper and safe operation of this device. Any other application, inconsistent with this manual, may lead to dangerous accidents. Restrict unauthorised personnel access to the device and train the personnel on handling the device. An operator of this device is any person who, by training, experience and knowledge of applicable reference standards, manuals and occupational health and safety regulations has been authorised for the essential operation with the device and who is capable of identifying and avoiding the hazards related to operation of this product.

Always append this Operating Manual with the device if transferred to a new owner. The Operating Manual contains detailed information about assembly, installation, commissioning, use, repairs and maintenance of the device. If the device is used as intended, this Manual will provide sufficient guidance to qualified personnel. Keep this Operating Manual close to the device and easily accessible at all times. As required by continuous improvement of the product, the manufacturer has the right to amend this Manual or make changes to the device without prior notice. Enbio Technology Sp. z o.o. shall not be liable for damage incurred during the wait for warranty service, any damage to the Customer's property other than this device, or errors caused by improper installation and/or improper operation of the device.

 Detailed recommendations, counter indications and warnings are described in the relevant sections.

## 2. SCOPE OF DELIVERY AND UNPACKING OF THE DEVICE.

### 2.1 Unpacking of the device

 If the steriliser was transported or stored at a temperature or humidity different than that at the location of installation, wait for 60 min. When moved from a cold room to a warm one, the device may contain moisture that, by negatively affecting the device's electrical components, may cause damage to it after startup.

 Remove the device from its packaging carefully.

 Attention! Check the packaging and its contents for external damage. If damage is found, contact the seller or the transport enterprise to prepare a damage report.

It is recommended to leave the box for possible autoclave transport.

### 2.2 Standard equipment

Verify the contents of the packaging in which the device has been delivered prior to installing it. The delivery packaging should contain:

1. ENBIO S/ ENBIO PRO steriliser
2. Water and condensate connection cables, rubber plugs for water/condensate containers
3. USB drive
4. Operating Manual
5. HEPA filter
6. TÜV Certificate
7. Validation report
8. ENBIO S/ ENBIO PRO certificate
9. Invoice

 According to Directive 2012/19/EU of the European Parliament and of


the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE) waste electrical and electronic equipment must not be disposed of or stored with household waste. Take the waste device to the nearest WEEE collection point registered with the General Inspectorate of Environmental Protection; this unit handles selective waste collection.

## 3. DEVICE INSTALLATION

 We recommend reading this user manual carefully before using the ENBIO S/ ENBIO PRO device. Follow all applicable safety guidelines and OHS regulations when operating the device.

 Mounting the HEPA filter. For reasons of transport safety, the HEPA filter has not been installed in the device. Remove it from the bag placed in the carton and tighten it yourself in a specially designated place on the back of the device. The filter should be screwed in manually until resistance is felt.

- a. The device should be positioned on a flat, even surface. Do not use the device if it is inclined.
- b. The device should be connected to a power supply that is earthed, equipped with fuses and has the same voltage rating as that indicated on the device.
- c. Demineralised or distilled water can be used in the device. Under no circumstances should tap water be used.
- d. Connect the connection tube included with the device to the water supply quick-release coupling on the device's rear panel, marked as WATER IN. Submerge the other end of the tube in the water supply container. The device is equipped with a water suction pump, there is no need to position the water container above or on the same level as the device. In order to secure the water supply tube, use the plug included in the delivery and place the plug in the opening of the water supply container. The minimum water load in the tank is 300 ml. Make sure that the blue hose of the inlet water is always immersed in the water.
- e. The wastewater formed after the water is turned into steam during the sterilisation process can be removed using the tube provided, which should be connected to the


port at the device's rear panel, marked as WATER OUT. The wastewater can be removed directly to the sewerage or to a special container intended for wastewater. If using container, place the tube end inside the container and secure the inlet with the plug provided. The tube must not be submerged in the wastewater.

f. The wastewater container or the sewerage drain must be located below the device.

g. If using wastewater containers, we recommend using containers of the same volume as those used for the deionised water. Emptying them concurrently with replacing/filling the deionised water containers will prevent overflow.

**⚠️ Correct positioning of cables in the water supply and wastewater containers.**

h. Leave 5 cm of space behind the device and 1 cm on each side from walls or other elements in order to ensure sufficient ventilation.

i. The device should be positioned in a way that ensures easy access to the main switch located on the rear panel of the device.

j. Do not position the device near washbasins or other places where it could be poured with water - possible short-circuit.

k. Install the device in a well-ventilated room, away from heat sources and rooms where mixtures of gases or liquids, and other hazardous agents may form.

l. Ensure the following environment conditions: operating temperature range +5°C to +40°C/ relative humidity 0–90%, storage temperature range from -20°C to +60°C/ relative humidity 0–90%.

**Enbio S and Enbio PRO devices are designed for self-assembly by the end user and do not require any special installation at the place of use. The user is responsible for the correct installation of the device on spot, according to this manual**

### 3.1 Water quality

ENBIO S/ ENBIO PRO sterilisers use demineralised or distilled water to form steam during the sterilisation process. The total mineral content in the water used for sterilisation must be lower than 10 ppm, or for conductance measurements, lower than 15 µS/cm.

Standard tap water has hardness within the 2–3 mmol/l range and must not exceed 5 mmol/l according to current regulations, making it unsuitable for use in ENBIO S/ ENBIO PRO sterilisers. Therefore tap water cannot be used as feed water for ENBIO S. / ENBIO PRO sterilizers.

The table below presents the hardness and conductance parameters of water used in steam sterilisation according to EN 13060.

Acceptable parameters of water used for sterilisation	
Hardness	< 0,02 mmol/l
Conductance (at 20°C)	< 15 µS/cm
Chemical additives	No chemical agents or additives must be added to the water used in the sterilisation process, even if they are intended specifically for use in steam generators, steam generation or for use as additives in sterilisation, disinfection, cleaning or corrosion protection.

**⚠️ Water conductance above 50 µS/cm may have a major impact on the sterilisation process and cause serious damage to the steriliser, and constitute ground for voiding the warranty. Use of water with impurities level exceeding the levels specified in the EN 13060 standard in the steam generator can significantly shorten the steriliser's lifetime.**

**⚠️ The distilled water in the supply tank should be replaced at least**

once every three months due to the increasing conductivity caused by prolonged contact with air. If the tank was contaminated, it should also be changed to a new one. The tank should be closed with the attached cap. Then the water does not change its properties so quickly.

 The manufacturer's warranty expires when the autoclave has been used with water services of a quality not compliant with the recommended one

## 4. TOOL PREPARATION AND LOADING

Only clean and dry tools may be sterilised. For this reason, before loading tools on the tray, clean and disinfect the tools in accordance with current regulations. Residue of agents used or solid particles may prevent the sterilisation process from completing successfully. Furthermore, sterilisation of tools not subjected earlier to pre-cleaning may cause damage to both the tools and the steriliser.

**During 134 C FAST program we do not use packages.**

If the tools were covered in grease, remove its excess.

Optimum method of positioning tools to be sterilised on the tray:

- For non-packaged tools – place the tools on the tray in such a manner so they do not contact each other directly. This will accelerate the drying process.
- For packed tools – position them on the tray in disposable sleeves as recommended by the pack manufacturer. Position packages with either the paper or film sides facing each other. Otherwise the packages may cement with each other during the sterilisation.

### 4.1. Tool pack preparation

#### 4.1.1 Characteristics of a sterilisation pack

It is recommended to use sterilisation packages that meet the requirements of the standards EN ISO 11607-1:2019, EN 868-2-10:2017-3.

A suitable pack is characterised by:

- good permeation of the sterilisation agent to the inside of the pack – resistance to damage during the sterilisation process,

- ensuring tight, durable sealing of the contents and their safe removal for further use,
- forming a barrier for microorganisms and undesired substances such as adhesive, ink from the label or a chemical test.

#### 4.1.2 Rules for arranging tools on a tray

Sterilized instruments should not protrude beyond the outline of the sterilization tray, special attention should be paid to sterilized instruments without packages. The tools must be positioned in such a way that no part of them falls into the holes of the tray, and does not rest on the edge of the sterilization tray or protrudes above the tray outline.

Failure to comply with the above recommendations may damage the sterilization chamber phase, which will result in the sterilizer being sealed.

- Sterilized instruments in packages: Arrange in a tray so that the package does not come into contact with the door seal and the phase of the sterilization chamber. Failure to comply may result in a lack of tightness in the device.
  - Do not exceed a maximum weight of 500 g for ENBIO S and 800 g for ENBIO PRO.
  - Special attention should be paid so that the ends of the packets do not protrude out of the sterilizer tray, which may cause the packet to jam during closing and lead to leakage of the sterilizer working chamber
  - It is recommended that when the working chamber is significantly loaded, the first packages should be directed with the foil side to the bottom of the tray.
- This guarantees faster and more efficient drying of packets.
- We do not use packages in the 134 C FAST program. Unpackaged sterilized instruments are intended for immediate use.

Not following the manufacturer's instructions will be associated with the loss of warranty on the device.

### 4.1.3 Principles of packing tools for sterilisation

#### Sterilisation pack type


Disposable paper and film packages

#### Principles of packing tools

- Bags should be filled only to 3/4 volume to allow proper sealing and minimize the risk of breakage
- a distance of 30 mm should be kept between the welding and sterilized equipment
- Protect sharp edges to avoid damage to the packaging
- the packaging material must not be laid loosely or stretched so that it does not affect pressure changes during sterilization
- the equipment should be stacked so that the paper side contacts the paper side as the sterilizing agent penetrates
- and air exchange can only take place through paper
- a label should be placed on the packaging with information about the content of the packaging, the code of the packer, date of sterilization and expiry date and sterilization parameters
- it is recommended to insert a sterilization strip into each process that discolours as a result of the correct sterilization cycle

Sample placing of sterilisation packages.

Example of solid batch


## 5. STARTING THE DEVICE

### 5.1 Commissioning

Before initiating the sterilisation cycle, turn the device on using the main switch located on the rear panel of the device. Make sure that water supply and wastewater cables are connected correctly, and that water is present in the water supply container, while the wastewater tank is empty, in order to prevent overflow. Monitor the water level in the tank regularly, depending on how frequently you perform your processes.

Place tools or materials in the working chamber tray and close the chamber and turn the knob locking the front of the device clockwise.

#### MAIN SWITCH

### 5.2 Program selection

Depending on the type of load to be sterilised, the user is responsible for selecting the appropriate program dedicated to the given type of load, in accordance with the manufacturer's recommendations for sterilisation.

PROGRAM ENBIO S	134°C FAST	134°C	121°C	PRION**
Type of load	Unpacked instruments	Packaged and non-packaged tools	Packaged and non-packaged tools	Packaged and non-packaged tools
Process temperature	134°C	134°C	121°C	134°C
Pre-vacuum number	1	3	3	3
Sterilisation duration	3,5 min	4 min	15 min	18 min
Drying duration	-	3 min 4 min <b>ENBIO PRO</b>	5 min 5 min <b>ENBIO PRO</b>	5 min
Total process duration*	100g: 7 min 100g: 10min <b>ENBIO PRO</b>	100g: 13 min 200g: 18 min <b>ENBIO PRO</b>	100g: 26 min 200g: 31 min <b>ENBIO PRO</b>	800g: 43 min
Max. water consumption	105 ml 140 ml <b>ENBIO PRO</b>	115 ml 190 ml <b>ENBIO PRO</b>	110 ml 180 ml <b>ENBIO PRO</b>	230 ml
Class	S	B	B	B

\*Ambient temperature can affect the process extension.

\*\*The duration of the first process may be longer due to the need for the device to heat up.


\*\*PRION program is available only in ENBIO PRO.

**The displayed time is approximate. The pre-vacuum time is updating and depends on the form and weight of the charge.**

Additionally, the Enbio PRO autoclave has a PRION program dedicated as one of the stages of decontamination of objects which are suspected to have had or may have had contact with pathologically altered prion proteins (e.g. Creutzfeldt-Jakob disease, BSE, etc.). Detailed information and recommendations on the control of transmissible spongiform encephalopathies are presented in the document "WHO Infection Control Guidelines for Transmissible Spongiform Encephalopathies - Report of a WHO consultation. (Geneva, Switzerland, 23-26 March 1999): It is the responsibility

of the user of the device to comply with the above guidelines.

**The 134°C program is recommended for the majority of sterilised materials due to the short duration of the entire program. The 121°C should be used to sterilise all other materials that cannot be subjected to sterilisation at the temperature of 134°C. Do not exceed a maximum weight of 500 g for ENBIO S and 800 g for ENBIO PRO.**

From here you can execute the Program, go to the Test, Information menu, and the COUNTERS menu.

In the Program menu, you can choose the 121°C temperature programs, 134°C FAST and 134°C PRION. When the chamber is opened, the **DOOR OPEN** symbol blinks.

When the chamber is locked by turning the lock knob clockwise, the display

will show the **START** symbol, indicating the chamber has been closed correctly.

Now you can select the program by pressing the symbol of the temperature you want to perform the

sterilisation in **121°C** or The selected program is initiated by pressing the corresponding symbol **START**

If the USB drive has not been inserted in the device, the USB symbol is not displayed in the bottom right corner of the screen, and a message about the missing USB drive is displayed.

The program data will not be saved. You can continue working without saving the data on the USB drive by

pressing the **YES** field or stop working by selecting the **NO** field to

insert the drive in the port and start the program from the beginning. If you decide to continue working or the **START** field has been selected, the screen will display a chart of pressure during the entire process, with the current stage of the program indicated, while information

on subsequent stages is displayed in the upper left corner of the screen. When a program is being run, the screen displays the temperature of the selected sterilisation program **121°C** or **134°C**, the current temperature of the process chamber in the bottom left corner **116.7°C**, the pressure currently in the chamber in the bottom right corner **0.30 Bar** while the process duration left is displayed in the upper right corner of the



screen **12 min left**. This is the expected time, which may be extended due to the mass and type of charge.

During the process, the **STOP** field is displayed instead of the **START** field, enabling you to stop the process at any time.

The upper left corner of the screen displays the status of individual subsequent program stages, e.g. –

- Chamber lock **Locking door.**
- Process chamber heating **Heating up**

**!** If the process was completed successfully, the display will alternately show information that the process has been completed and the load is sterile, and that the chamber may be opened. After the process in 134C FAST program tools are wet and hot.


By pressing the **FINISH** field, you return to the start screen.

**!** **ATTENTION!** When the process is completed, the chamber, the tray and the load are hot. Maintain particular care and use protective gloves to remove the load, or wait until it cools.

Performing the sterilisation process in the ENBIO S/ ENBIO PRO device does not affect material biocompatibility. All components of the device that are in direct contact with the sterilised load have no toxic, sensitising or irritating effects.

### 5.3 Test programs


By pressing the **Test** field, you can go to the test program menu.

From here you can select the vacuum leak test program and the Helix/B&D test program. Select the appropriate program by pressing the relevant field on the display.

When the process chamber is closed, the **DOOR OPEN** information changes to **START** and by pressing this field you launch the selected test program.

If the USB drive has not been inserted in the device, the USB symbol is not displayed in the bottom right corner of the screen, and a message about the missing USB drive is displayed.

The test program data will not be saved. You can continue working without saving the data on the USB drive by pressing the **YES** field or stop working by selecting the **NO** field to insert the drive in the port and start the test program from the beginning.

Program ENBIO S / ENBIO PRO	Bowie & Dick / Helix	Vacuum leak test
Process temperature	134°C	–
Pre-vacuum number	3	1
Sterilisation duration	3,5 min	–
Drying duration	3 min	–
Total process duration	15 min	16 min


### Vacuum leak test

The vacuum leak test may only be performed on a cold device, before work is commenced. The vacuum leak test enables testing the autoclave for the presence of leaks. The following are checked during the test:


- Vacuum pump performance.
- Pneumatic system sealing.

When the vacuum leak test program is selected and launched with the **START** button, the vacuum leak test progress screen is displayed.


Information on pressure loss in the process chamber, and the test duration are displayed.

When the test program is completed, the following screens are displayed alternately.

When the test program was completed successfully.

When the test program was not completed successfully.

After pressing the  field, the start screen is displayed.

 The process chamber must be completely dry and cold during the vacuum leak test. Otherwise, the vacuum leak test results may not be fully reliable, even if the steriliser is fully functional. When the test is completed, a message with the results will be displayed. If the result is negative, check, clean or replace the seal, clean the front edge of the chamber and repeat the test.

If the device fails the test again, contact your supplier or the manufacturer.

## Bowie&Dick test

Perform the Bowie&Dick test daily before commencing work, in order to verify that the device performs sterilisation correctly.

The Bowie&Dick test, also known as the steam penetration test, imitates a small, highly porous load.

It contains sheets of paper packaged in a small pack containing a chemical indicator (a physicochemical test).


This test assesses the device's performance in sterilising charges composed of porous objects:

- Pre-vacuum performance and steam penetration.
- Saturated steam temperature and pressure, reached for a specific time.

How to perform the test:




- Perform the test with an empty chamber, as per the EN 13060 standard.
- Place the Bowie-Dick test pack in the chamber, in the middle of the tray.

When the Helix/B&D test program is selected and launched with the  button, the program progress screen is displayed.

Process parameter information is displayed.

The Helix/B&D test program can be stopped at any time by pressing the  field.

When the test program is completed, the following screens are displayed alternately.

You can safely open the steriliser's process chamber.

When the process chamber is opened, the start screen is displayed.

- Remove the test pack.

 **WARNING! The package will be hot.**

In order to interpret the test correctly, read the instruction provided by the test pack manufacturer.

- Open the pack and remove the chemical indicator from inside.

**POSITIVE RESULT**  
chemical indicator has changed color on a dark uniform over the entire surface.

#### NEGATIVE RESULT

In the middle of the test field remains clear because of the residual air in the middle of the device under test.

Any change of colour, uneven colouring of the test, indicate the presence of air during the test cycle, caused by faulty operation of the steriliser.

#### Helix test

The Helix test represents sterilisation of tools with A-type holes in accordance with the EN 13060 standard. It consists of a 1500 mm long tube open on one side, and a closed test capsule on the other. The indicator strip is located in the test capsule.

Helix test set

This test is used to assess the device's performance in sterilising recessed and porous charges, in particular:

- Pre-vacuum performance and the rate and uniformity of steam penetration.
- Saturated steam temperature and pressure, reached for a specific time.

How to perform the test:

- Perform the test with an empty chamber, as per the EN 13060 standard.
- Place the test strip in the capsule. Read the instruction provided by the test's manufacturer.
- Close the capsule.
- Place the test in the middle of the tray in the chamber.
- After completing the test, open the steriliser and remove the test.


 **WARNING! The package will be hot.**

In order to interpret the test correctly, read the instruction provided by the test pack manufacturer.

– Open the capsule and remove the test strip.

**POSITIVE RESULT**  
All fields of the indicator strip became dark.

**NEGATIVE RESULT**  
Part of the indicator strip did not change color on the dark due to the presence of air inside the capsule. Part of the indicator strip did not change color on the dark due to the presence of air inside the capsule.

A part of the indicator strip did not change colour to the dark one due to the presence of air in the capsule.

An insufficient change of colour of the indicator strip indicates the presence of air during the test cycle, caused by faulty operation of the steriliser.

## 5.4 Information menu

The information menu can be accessed by pressing the  field.

Here, information about the device type, serial number, number of performed processes, amount of free memory available on the USB drive for saving process data, and the service menu  with the process counters for sealing, filter and the next service inspection can be displayed. You can also change the date and time.

In order to set the date or time, touch the digits on the display. When a field is selected, it starts to blink  and arrows used to change values are displayed, up  and down . This way you can correctly set the date and time.

Pressing the number again confirms it, and you can change another parameter.

We can choose the language in the same way by clicking on the shortcut. The button marked with the letter B turns off and on the blue backlight in the back of the screen

Clicking the LED button launches the backlight control menu, which is located on the sides of the device. LED lighting has two modes:

- Free mode, in which the user (by moving the sliders) sets the colors, intensity and brightness of the light.
- Continuous mode, which indicates the stages of the entire sterilization process with colors

**LED lighting is available only in ENBIO S devices with colored housings: Midnight Blue, Dolly Pink, Ashy Stone, Yellow Sunflower.**



### 5.4.1 Counters

The ENBIO S/ ENBIO PRO steriliser counts the number of performed processes and uses it to notify you about the recommended dates of replacing elements subject to wearing down, and about required service inspections.

No	Name	Recommended frequency of replacement (cycles)	Yellow (Nearing the replacement date, number of cycles)	Red (Exceeded replacement date, number of cycles)
1.	HEPA filter	after 1000	from 980	after 1000

By pressing the **COUNTERS** field, you go to the counter display.

The number of performed processes is on the left hand side, and on the right hand side – the number at which the given element should be replaced or a service inspection performed. **980/1000**. After replacing a filter or seal, the values can be reset by the user by pressing the **RESET** button. The service inspection value can only be reset by an authorised service.

When nearing a value when replacement of an element or a service inspection is recommended, the values will be highlighted in yellow.

If the limits are exceeded, the value will be displayed in red.

During regular operation, info screens concerning replacement of individual elements or required service inspection are displayed alternately.

Counter values displayed in yellow or red do not prevent the device from operation. However, exceeding the required inspection date may significantly affect the device's operation and the load sterilisation process.


For replacing individual elements, please contact the manufacturer or supplier.

### 5.5. Restarting

Restarting the device is forced when a process is aborted by the user by pressing the **STOP** field, or when power or water supply is lost.

STOP

If the **STOP** field is selected, the following messages will be alternately displayed, notifying you that the process has been aborted by the user and pressure is being equalised in the process chamber, and a message notifying you that the process has not been completed correctly and the load is not sterile.




When the pressure is equalised in the process chamber, the following messages will be displayed alternately on the screen. You can freely open the device now.

The following screen will be displayed when the chamber is opened.

By selecting the **RESTART** field, you can return to the start screen.

However, in order to do so, you must enter the 4-digit security code **0000**.

If the code is entered incorrectly, a message will be displayed on the screen.

Enter the code again. The arrow enables cancelling incorrectly entered digits.

When the code is entered correctly, the start screen will be displayed.

## 6. MAINTENANCE AND CARE

### Tray cleaning

Maintaining tray cleanliness aids in maintaining correct functioning of

the device.

It is recommended to clean the internal part of the tray once a week using a mild, chlorine-free detergent that does not react with aluminium. After


cleaning, the tray must be thoroughly washed with water.

Dry the tray before reinstalling the tray and push it over the front face pins and push it down gently to lock it.

### **Cleaning the process chamber**

Maintaining the chamber cleanliness aids in maintaining correct functioning of the device.

It is recommended to clean the process chamber interior once a week using a mild, chlorine-free detergent. After cleaning, wipe the chamber with a soft cloth until dry.

To clean the tray well it must be removed from the front of the device. To do this, lift the tray gently up and pull it away from the front. The mounting studs have notches in which the drawer fits.

### **Cleaning external surfaces**

The external parts of the device should be cleaned using a soft cloth slightly wet with water and a mild detergent (chlorine-free and not reacting with plastics, varnishing coats and aluminium). Do not use strong detergents.

Use of mild detergents for maintaining the device does not affect the possibility of hazard related to toxic agents forming in contact with elements of the device.

### **Cleaning the seal**

It is recommended to clean the seal after 100 performed processes. Use warm water and a microfibre cloth (microfibre with silver particles is acceptable) to clean the device. Use of dull or sharp cleaning tools is not acceptable. Cleaning with chemical agents is not acceptable. Perform the cleaning when the device has cooled down, after opening the drawer. Maintain caution and do not bend the drawer. After cleaning, leave the device open until the seal dries. During this time, protect the device from damage. After cleaning and drying the seal, it can be lubricated with a silicone lubricant.

### **Replacement of elements subject to wearing down**

Elements subject to wearing down should be replaced periodically to ensure failure-free operation of the steriliser.

A message on the screen will notify the user when individual elements should be replaced.

During regular operation, info screens concerning replacement of individual elements or required service inspection are displayed alternately. They are described in detail in the „Warning messages and error codes“ section.

## Cleaning the water container

In order to ensure correct parameters of the water supplying the device, it is recommended to check the water tank at least once a quarter. If contamination is found, the tank should be drained, cleaned and refilled with new water.

 In order to ensure efficient sterilisation and correct operation of the device, it is recommended to observe the replacement dates for elements subject to wearing down.

### 6.1 Replacement parts

The following table includes elements subject to periodic replacement, and elements subject to natural wear and tear. Replacement parts should be ordered directly from the manufacturer. Use of other replacement parts voids the warranty and does not guarantee correct functioning of the device.

Name	Part no.
Front seal	ST1-UL1
Bacteriological filter	DZ0035
Connection/water supply tube	ST1-HW1
Connection/water supply tube	ST1-HW2
Rubber plug for the water container	ST1-KS1
Rubber plug for the condensate container	ST1-KS2

### 6.2 Periodic inspections

In order to ensure correct functioning of the ENBIO S/ ENBIO PRO steriliser, it is recommended to perform periodic service inspections and replace parts subject to wearing down in accordance with the following schedule, and periodic inspection of individual steriliser elements in accordance with the following guidelines.


Name	Replacement frequency
Bacteriological filter	every 1000 cycles or every 12 months
Connection/water removal tube	if damage is observed
Plugs for water/condensate containers	if damage is observed

Element subject to inspection	Inspection frequency
Front seal	weekly or in the event of incorrect functioning – performed by the user
Bacteriological filter	weekly – performed by the user
Connection/water removal tube	weekly or in the event of incorrect functioning – performed by the user
Container plugs	weekly – performed by the user

## 7. DATA ARCHIVING

The progress of each performed sterilisation is automatically saved on a data carrier (USB drive). The data can be used only for archiving, the sterilisation process correctness is directly communicated by the device.

The USB port is located on the rear panel of the device. It is recommended to periodically archive the data on another carrier, e.g. a desktop PC, laptop.

## 8. ENBIODATAVIEWER

The EnbioDataViewer software enables viewing and archiving sterilisation programs on a computer, and printing them.

Minimum requirements for installing the software:

Operating system – Windows – at least Windows 7 or newer

Free disk space – at least 100 MB

Minimum CPU parameters – at least 1 GHz

Minimum RAM – at least 512 MB


Screen resolution – at least 1200 x 720 or more


The software is delivered with the device and can be found on a removable disk - the pendrive or the latest version can be downloaded from the manufacturer's website <http://enbio-group.com/steamjet/steamjet-software-serwis/>

## 8.1 Software installation


To install the software, double click on the software installation file. After performing this operation, the installation window regarding the language selection will be displayed


After confirmation, you must accept the license terms for the installed software.


Next, the information about placing the software shortcut on the computer desktop will be displayed.

After making your selection, click „Next“. By clicking the Install button you will install the EnbioDataViewer software.


After installation, the following message is displayed.


We can now run the software or finish the installation without running the software by clicking the Finish button.

The main program window is displayed.

## 8.2. Program construction and main functionalities

The main window consists of three main areas

All processes, which have been synchronized from pendrive were sorted after the dates of performance

Temperature and pressure diagram together with the main data regarding autoclave and process (date of completion).

Data on the subsequent stages of the process.

the most important parameters of sterilization.

The ability to save a note for each process.

The dark blue color has been marked with function keys, e.g. „PDF Report” that will allow you to print a protocol from the process.



Drop-down menu:

By clicking on the File window we have access to the options:

- Loading the saved process flow from the memory of the pendrive or from another location
- Printing a saved program
- Implementation of the report to a PDF file
- Exporting data to the database to get it sent to the producer in case of problems
- Exporting data to CSV format
- Closing the program


By clicking on the Tools window we have access to the options:

- Synchronization of all files with saved programs from the memory of the pendrive
- Search for any saved process from the database
- Adding your own logo to PDF reports

By clicking on the Help drop-down menu, you have access to the options:


- About program

## Search


The program allows you to search for processes by:

- Range of dates
- Sterilization number
- The type of process
- Result

## PDF report

The program allows you to generate a report from every process performed by the autoclave. It contains all necessary process data and the result of sterilization.


## 9. WARNING MESSAGES AND ERROR CODES

If any irregularities in the device's operation occur, the screen will display relevant warning messages and error codes.

### 9.1 Warning messages

The element subject to replacement is highlighted in red, the screens are displayed alternately.

Screens concerning seal replacement, with the number of processes remaining until replacement.

Screens concerning filter replacement.

### 9.2 Information codes

Screen for overpressure or underpressure resulting from natural chamber cooling processes.

Message resulting from interruption of the process after the sterilization stage - during drying.

### 9.3 Error codes

The following table includes the error codes that may be displayed during the use of the steriliser.


Error code	Description	Recommendations
1	"Chamber over temperature"	Maximum temperature in the chamber exceeded Contact the service
2	"Steam gen. over temperature"	Excessive steam generator temperature Contact the service
3	"Process over temperature"	Excessive process temperature Contact the service
4	"Overpressure error"	Pressure error Contact the service
5	"Sterilization pressure too low"	Pressure too low during sterilisation Check water level and connection. Contact the service
6	"Sterilization temp. too low"	Sterilization temperature too low Check water level and connection. Contact the service
7	"Too high pressure during drying"	Pressure too high during drying Check if the outlet tube is not submerged in water. Contact the service
8	"Too many steam pulses/no water"	Too many steam impulses. No water supply. Check water supply level and tubing connections. Contact the service
9	"Drainage error"	Drainage clogged Check wastewater level and tubing connections. Contact the service
10	"Chamber heating error"	Chamber heating error Contact the service
11	"Steam generator heating error"	Steam generator error Contact the service
12	"Prevacuum fail/check condensate outlet"	Vacuum pump/drainage error Check wastewater level and tubing connections. Contact the service

13	"Power failure"	Temporary power loss during operation Confirm error.	
14	"Pressure during standby"	Pressure exceeded during standby Confirm error. Contact the service	
15	"Locking door error"	Door lock error Contact the service	
16	"Unlocking door error"	Door unlocking error Contact the service	
17	„Valve V3 / HEPA filter error“	V3 valve / HEPA filter error Check filter cleanliness/replace filter. Contact the service	
18	"Pressure sensor error"	Pressure sensor error Contact the service	
19	"USB disc error / Change disc"	Write error on the pendrive - media damage Write error on the pendrive - damage to the media. Rip content from the current pendrive - purchase of a new one and use of a new one	
31	"Internal flash error"	Memory error Contact the service	
<b>Info</b>			
	„Aborted by user“	Process aborted by the user. Non-sterile insert if interrupted during or before the sterilization process.	
	„Vacuum test failed“	Vacuum leak test error Contact the service	
	„No USB memory“	No USB drive Check the USB port, insert the drive. Contact the service	
	„Equalizing pressure“	Pressure during stoppage. Equal to atmospheric pressure The message occurs in specific cases as a result of natural processes. In the case of a frequent appearance of a message, contact the service.	

Sample error codes are presented below.

Process aborted by the user.

Screens displayed alternately: equalising pressure, please stand by.

The process has not been completed correctly. The load is not sterile.

Vacuum leak program error. Error screen: you can continue working.

Error no. 5  
Pressure in the process chamber too low.

The process has not been completed. The load is not sterile. Equalising pressure in the process chamber.

## 10. WARRANTY CLAIM HANDLING

In order to report a problem with the device, fill out the Warranty Claim Form available on the manufacturer's website [www.enbio.com](http://www.enbio.com). Our Technical Service will contact you as soon as possible. If the device has been damaged during transport, file your warranty claim with the delivery note and photographic evidence of the damage found.

**⚠ ATTENTION!** The warranty claim process will begin once our Technical Service receives a properly completed Warranty Claim Form.

If you send the device to the Technical Service, clean the chamber and tray, perform decontamination and correctly secure the device for transport. Preferably send the device in the original packaging. If you lack an appropriate packaging, please contact the Technical Service or your supplier.

If the device needs to be transported:

- Disconnect the demineralised water and condensate tubing.
- Wait until the process chamber cools down.
- Use the original or other appropriate packaging with protection lining.

The sender bears complete liability for damage during transport to the Technical Service.

## 11. WARRANTY TERMS AND CONDITIONS

ENBIO sterilisers are covered by a standard 12-month warranty. Detailed warranty terms and conditions are available from the supplier of this device.

## 12. TECHNICAL INFORMATION

Technical data	ENBIO S	ENBIO PRO
Power supply	230 V/50Hz	230 V/50Hz
Installed power	3,25 kW max	3,25 kW max
Operated Power	1,6 kW	1,6 kW
Maximum electric current consumption	15 A	15 A
Operating pressure	2,1 bar	2,1 bar
Maximum pressure	2,3 bar max	2,3 bar max
Maximum process temperature	137°C	138°C
Process chamber capacity	2,7 l	5,3 l
Mass	15 kg	20 kg
Process chamber dimensions (LxWxH)	292 x 192 x 45 mm	300 x 200 x 90 mm
External device dimensions (LxWxH)	561 x 252 x 162 mm	561 x 270 x 202 mm
Protection rating	IP20	IP20
Noise level	49dB(A)	49dB(A)
Process data archiving	USB drive	USB drive


### Surrounding conditions:

Operating temperature range	from +5°C do +25°C
Relative humidity	0-90%
Storage temperature range	from -20°C do +60°C
Relative humidity	0-90%
Surrounding pressure range	900-1100 hPa

Rating plate located on the bottom of the device.

enbio		CE 2274		RoHS	
REF	Enbio S	TYPE	B	SN	ST01-PL-20-03312
<b>Power supply</b> 220-240 V AC 10A 50/60Hz 2,25 kW max.			<b>Sterilization chamber</b> Max. pressure 2,1 bar Min. pressure -0,90 bar Max. temperature 134 °C Pressure test 9,4 bar Chamber volume 2,7 dm <sup>3</sup>		
<b>Manufacturer</b> Enbio Technology Sp. z o.o. ul. Sportowa 8 C, 81-300 Gdynia, POLAND tel.: +48 22 469 81 99 www.enbio.com			<b>Maxima</b> Max. working pressure 2,3 bar Max. working temp. 137 °C		

enbio		CE 2274		RoHS	
REF	ENBIO PRO	TYPE	B	SN	ST02 - PL - 19 - XXXXX
<b>Power supply</b> 220-240 V AC 15A 50/60Hz 3,25 kW max.			<b>Sterilization chamber</b> Max. pressure 2,1 bar Min. pressure -0,90 bar Max. temperature 134 °C Pressure test 9,4 bar Chamber volume 5,3 dm <sup>3</sup>		
<b>Manufacturer</b> Enbio Technology Sp. z o.o. ul. Sportowa 8 C, 81-300 Gdynia, POLAND T. +48 22 469 81 99 www.enbio.com			<b>Maxima</b> Max. working pressure 2,3 bar Max. working temp. 137 °C		

Test connector - to be used by authorized service only. If it is determined that the user has used it, the warranty will be voided.

## 13. EC DECLARATION OF CONFORMITY

### Declaration of Conformity RoHS

Company: Enbio Technology Sp. z o. o., ul. Sportowa 8 C, 81-300 Gdynia, Poland

declares with sole responsibility, that medical devices: **STEAM STERILIZERS**, models:

- **Enbio S**
- **Enbio PRO**
- **Enbio XS**

complies with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

Gdynia  
14.01.2021

Marek Krajczyński

A handwritten signature in blue ink, appearing to read "M. Krajczyński".

President of the Board

Sebastian Magrian

A handwritten signature in blue ink, appearing to read "S. Magrian".

Member of the Board

**Dane producenta**

**Enbio Technology Sp. z o.o.**

ul. Sportowa 8 C

81-300 Gdynia, Polska

T. +48 22 469 81 99

[www.enbio.com](http://www.enbio.com)