

GLASS FUNNEL



DESCRIPTION

The Teqler glass funnels are used for filling and transferring liquids from bottles and vessels with small openings. The borosilicate glass used in production guarantees high chemical and temperature tolerance up to 140°C. You can also use filter paper with this funnel.

PRODUCT DETAILS

- Teqler glass funnel
- **Made of borosilicate glass**, clear
- For transferring liquids and fine substances
- **Suited for filtration (with filter paper)**
- High chemical and temperature tolerance (up to +140°C)
- 3 different sizes available
- 1 piece

Ø Fill Opening	Spout Length	Ø Spout Opening
----------------	--------------	-----------------

40 mm	40 mm	ca. 5 mm
60 mm	60 mm	ca. 7 mm
100 mm	90 mm	ca. 10 mm

Overview


Article Number:
T135851-T135853

Category:
Lab Equipment



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
 8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09
  +352 267149 19
  info@netmed.lu

Social Media

