



germanengineering

PO 30 Pulsoximeter

Ermittlung der arteriellen Sauerstoffsättigung (SpO₂) und der Herzfrequenz (Pulsfrequenz)

Sehr einfache und vollkommen schmerzfreie Messung

Klein und leicht für zu Hause und unterwegs

Besonders geeignet für Personen mit:

- Herzinsuffizienz
- Chronisch obstruktiven Lungenerkrankungen
- Asthma bronchiale

Anwendbar bei Sport in großen Höhen (z.B. Bergsteigen, Skifahren und Sportfliegerei)

Leicht ablesbares Farbdisplay

Display mit 6 Ansichtsperspektiven

Einstellbare Display-Helligkeit

Graphische Pulsanzeige

Batteriewechselanzeige

Abschaltautomatik

Inkl. Halteband und Gürtel-Tasche

Medizinprodukt

Inkl. 2 x AAA Batterie

Produktmaße: 61 x 36 x 32 mm

Produktgewicht: ca. 57g

3 Jahre Garantie

VE: 4 / Transportkarton: 8 x 4

EAN-Nr.: 421112545430 2

Art.-Nr.: 454.30

PO 30 Pulse Oximeter

Determining arterial oxygen saturation (SpO₂) and heart rate (pulse frequency)

Simple and completely pain-free measurement

Small and light – for use at home and on the move

Particularly suitable for persons with:

- Heart failure
- Chronic obstructive pulmonary diseases
- Bronchial asthma

Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)

Easy-to-read colour display

Display with 6 available views

Adjustable display brightness

Graphical display of pulse

Low battery indicator

Automatic switch-off

Incl. retaining strap and belt pouch

Medical product

Incl. 2 x AAA batteries

Product measurements: 61 x 36 x 32 mm

Product weight: approx. 57g

3 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 421112545430 2

Item no.: 454.30