



germanengineering



## PO 80 Pulsoximeter

Ermittlung der arteriellen Sauerstoffsättigung (SpO<sub>2</sub>) und der Herzfrequenz (Pulsfrequenz)

Sehr einfache und vollkommen schmerzfreie Messung

Klein und leicht für zu Hause und unterwegs

Besonders geeignet für Personen mit:

- Herzinsuffizienz
- Chronisch obstruktiven Lungenerkrankungen
- Asthma bronchiale
- Schlafapnoe

Anwendbar bei Sport in großen Höhen (z.B. Bergsteigen, Skifahren und Sportfliegerei)

Alarmfunktion - Grenzwerte individuell einstellbar

Kontinuierliche Aufzeichnung bis zu 24h

PC-Software zur detaillierten Darstellung und Auswertung  
Betriebssysteme: - Win XP - Win Vista - Win 7

USB-Kabel und USB-Ladegerät

Aufladbar per USB (integrierter Lithium-Ionen-Akku)

XL Farbdisplay mit 4 Ansichtsperspektiven

Einstellbare Display-Helligkeit

Graphische Pulsanzeige und zuschaltbarer Pulston

Batteriestandsanzeige

Abschaltautomatik

Inkl. Halteband und Gürtel-Tasche

Medizinprodukt

Produktmaße: 57 x 32 x 30 mm

Produktgewicht: ca. 42 g

3 Jahre Garantie

VE: 4 / Transportkarton: 8 x 4

EAN-Nr.: 421112545440 1

Art.-Nr.: 454.40

## PO 80 Pulse Oximeter

Determining arterial oxygen saturation (SpO<sub>2</sub>) and heart rate (pulse frequency)

Simple and completely pain-free measurement

Small and light – for use at home and on the move

Particularly suitable for persons with:

- Heart failure
- Chronic obstructive pulmonary diseases
- Bronchial asthma
- Sleep apnoea

Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)

Alarm function - limit values can be individually set

Continuous recording of up to 24h

PC software for detailed display and evaluation  
Operating system: - Win XP - Win Vista - Win 7

USB cable and USB charger

Can be charged by USB (integrated lithium-ion battery)

XL colour display with 4 available views

Adjustable display brightness

Graphical display of pulse and optional pulse tone

Battery status indicator

Automatic switch-off

Incl. retaining strap and belt pouch

Medical product

Product measurements: 57 x 32 x 30 mm

Product weight: approx. 42 g

3 year guarantee

Sales unit: 4 / Shipping carton: 8 x 4

EAN no.: 421112545440 1

Item no.: 454.40