

## Specifications : BT-410F Fixed Beamlight Headlamp

<b>Functional Characteristics</b>	
<b>LED Light</b>	
Illumination (at 25cm)	15,000 ~ 30,000lx
Color temperature	6,000Kelvin
Life time	50,000hrs
<b>Function</b>	
Battery state display	Green / red(low battery)
Battery module attachable	
Light Intensity Adjustable	
<b>Power</b>	
Adaptor	Input : AC 100 ~ 240V (50/60Hz)
	Output : DC 4.2V/750mA
Rechargeable bBattery	3.7V, 2,600mA Li-ion
	Operating time : 4hrs
	Charging time : 4hrs
<b>Standard Configurations</b>	
Operation manual	1ea
Battery module	1ea
Adapter	1ea
<b>Optiops</b>	
Loupe	x 3.5 (diopter)
Battery module	
<b>Physical Characteristics</b>	
<b>Dimension</b>	
Headband	534 ~ 638mm
Packing	282(W) x 187(D) x 117(H)mm
Carton box (5ea)	600(W) x 300(D) x 215(H)mm
<b>Weight</b>	
Headband	119g
Battery module	69g
Packing	670g
Carton box (5ea)	4.35Kg
<b>Environmental Conditions</b>	
Operating temperature	10 ~ 40°C (50 ~ 104°F)
Operating humidity	5 ~ 85% non-condensing
Storage temperature	-20 ~ 60°C (-4 ~ 140°F)
Storage humidity	0 ~ 95% non-condensing
<b>Warranty</b>	
1year	
<b>Certificates</b>	
KFDA, CE, FDA	

