

## SHARPS BINS | 48 LITER



### DESCRIPTION

The XL safety container from Teqler is pierce-proof, stable and equipped with a large opening that can be temporarily or irreversibly closed. Therefore, the container is ideal for discarding sharp and pointed objects, such as cannulas, complete syringes and disposable instruments. Thanks to a transparent lid, the filling level of the container with a volume of 48 liters can be checked at any time.

### PRODUCT DETAILS

- Discard container for sharp and pointed medical consumables
- Particularly puncture-resistant and stable
- Can be closed temporarily or irreversibly
- Extra large opening
- Transparent lid allows control of the filling level
- Filling volume: 48 liters
- Color: yellow
- Dimensions: 41 x 28,5 x 57 cm
- 1 piece

### Overview

Article Number:  
**T140348**




Category:  
**Sharps bins Collection**

Dimensions:  
**41 x 28,5 x 57 cm**



Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.  
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  [info@netmed.lu](mailto:info@netmed.lu)

Social Media

