



Typhoid Antigen Rapid Test Kit (Colloidal Gold Method)

The Typhoid Antigen Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection of Salmonella Typhi (*S. typhi*) antigens in human fecal, serum or plasma. This method is easy to operate, and the result can be obtained by observing whether there is a red and purple line in 15-20 minutes. It is intended to be used as a screening test and provides a preliminary test result to aid in the diagnosis of infection with *S. typhi*.

Test Procedure

- **Specimen:** Fecal, serum or plasma.
- **Specimen collection:** Specimens of all types are collected by conventional methods.
- **Test method:**



Specifications

Components	Test card/ Sample extraction tubes /Sample diluent/Sample droppers
Pack Size	25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C~30°C ;Valid for 24 months
Test samples	Fecal, Serum & Plasma
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Droppers	Sample Type
Typhoid Antigen Rapid Test Kit (Colloidal Gold Method)	BSK15S1S	25 Tests/Kit	1 × 3 mL/Tube	25	Fecal, Serum & Plasma
	BSK15M1S	50 Tests/Kit	2 × 3 mL/Tube	30	