


SURGICRYL[®]

MONOFILAMENT

POLYDIOXANONE



Type	Monofilament
Composition	Polydioxanone, a polymer made from polyester poly (p-dioxanone)
Coating	None
Colour	Violet
Tissue reaction	Minimal
Absorption	The hydrolytic action by which the material is broken down results in total absorption in approximately 180 to 210 days. Approximately 50% of tensile strength remain after 42 days.
Presentation	Box with 12 sutures



*For quality
&
safety*

www.sutures.be

SURGICRYL[®]

MONOFILAMENT

POLYDIOXANONE

CHARACTERISTICS

- High and long-lasting tensile strength
- Synthetic
- Breaks down by hydrolysis
- Monofilament with microscopic uniform diameter
- Pliability and softness
- Hermetically sealed packing

INDICATIONS

General soft tissue closing especially in

- General surgery
- Gastrointestinal surgery
- Ophthalmic surgery
- Orthopaedics
- Abdominal wall closure
- Gynaecology
- Plastic surgery
- Urology

BENEFITS

Outstanding strength retention for extended wound support

Minimal tissue reaction

Predictable absorption rate

Soft passage through the tissues without sawing, tissue drag and trauma

Absence of capillarity (no support of bacterial growth)

Excellent knot tie down

Good handling properties

Easy and secure knotting

Guaranteed seal and product sterility

CONTRAINDICATIONS

Cardiovascular surgery


Neurosurgery

Prosthetic devices


Being absorbable not to use when wound support for more than 6 weeks is required.

Caution in patients with delayed wound healing (e.g. compromised immune system).

TENSILE STRENGTH



ABSORPTION



 Without needle	Metric	USP	Length	Order number
	5	2	150 cm	<i>violet</i> 13500015

 Without needle	Metric	USP	Length	Order number
--	--------	-----	--------	--------------

Ligapack















2	3/0	250 cm	<i>violet</i> 13200025
---	-----	--------	------------------------






 Round bodied taper point	Metric	USP	Length	Order number
--	--------	-----	--------	--------------




1/2 circle




















HR 13 mm		0.7	6/0	75 cm	<i>violet</i> 13070113
		1	5/0	75 cm	<i>violet</i> 13100113
HR 17 mm		0.7	6/0	75 cm	<i>violet</i> 13070117
		1	5/0	75 cm	<i>violet</i> 13100117
		1.5	4/0	75 cm	<i>violet</i> 13150117
		2	3/0	75 cm	<i>violet</i> 13200117
HR 22 mm		1.5	4/0	75 cm	<i>violet</i> 13150122
		2	3/0	75 cm	<i>violet</i> 13200122
HR 26 mm		1.5	4/0	75 cm	<i>violet</i> 13150126
		2	3/0	75 cm	<i>violet</i> 13200126
		3	2/0	75 cm	<i>violet</i> 13300126
HR 30 mm		3.5	0	75 cm	<i>violet</i> 13350126
		2	3/0	75 cm	<i>violet</i> 13200130
HR 36 mm		3	2/0	75 cm	<i>violet</i> 13300136
		3.5	0	90 cm	<i>violet</i> 13350136
HR 40 mm		4	1	90 cm	<i>violet</i> 13400136
		3	2/0	75 cm	<i>violet</i> 13300140
		3.5	0	90 cm	<i>violet</i> 13350140
HR 48 mm		4	1	90 cm	<i>violet</i> 13400140
		5	2	90 cm	<i>violet</i> 13500140
		3.5	0	90 cm	<i>violet</i> 13350148
		4	1	90 cm	<i>violet</i> 13400148
		5	2	90 cm	<i>violet</i> 13500148




 Round bodied taper point	Metric	USP	Length	Order number	
	DR 13 mm 	1	5/0	75 cm	<i>violet</i> 13100213
		1.5	4/0	75 cm	<i>violet</i> 13150213
	DR 17 mm 	1	5/0	75 cm	<i>violet</i> 13100217




 Round bodied taper point double armed	Metric	USP	Length	Order number	
	2 x HR 13 mm 	1	5/0	75 cm	<i>violet</i> 13100113-2
	2 x HR 17 mm 	1	5/0	90 cm	<i>violet</i> 13100117-2
		1.5	4/0	75 cm	<i>violet</i> 13150117-2
	2 x HR 26 mm 	2	3/0	90 cm	<i>violet</i> 13200126-2

 Round bodied taper point double armed	Metric	USP	Length	Order number	
	2 x DR 13 mm 	0.7	6/0	75 cm	<i>violet</i> 13070213-2

 Round bodied taper point loop	Metric	USP	Length	Order number	
	HR 40 mm 	3.5	0	150 cm	<i>violet</i> 13350940
		4	1	150 cm	<i>violet</i> 13400940
HR 50 mm 	3.5	0	150 cm	<i>violet</i> 13350950	
	4	1	150 cm	<i>violet</i> 13400950	
	5	2	150 cm	<i>violet</i> 13500950	
HR 65 mm 	4	1	150 cm	<i>violet</i> 13400965	

 Reverse cutting		Metric	USP	Length		Order number
	DS 12 mm 	0.7	6/0	45 cm	<i>violet</i>	13071512
		1	5/0	45 cm	<i>violet</i>	13101512
		1.5	4/0	75 cm	<i>violet</i>	13151512
	DS 16 mm 	1	5/0	75 cm	<i>violet</i>	13101516
		1.5	4/0	75 cm	<i>violet</i>	13151516
	DS 16 mm Easy Pass 	1	5/0	50 cm	<i>violet</i>	13101516P
	DS 16 mm black 	1	5/0	75 cm	<i>violet</i>	13101516VB
	DS 19 mm 	1	5/0	45 cm	<i>violet</i>	13101519
		1	5/0	75 cm	<i>violet</i>	13101519-75
		1.5	4/0	75 cm	<i>violet</i>	13151519
		1.5	4/0	90 cm	<i>violet</i>	13151519-90
	DS 19 mm black 	2	3/0	75 cm	<i>violet</i>	13201519
		1.5	4/0	75 cm	<i>violet</i>	13151519VB
		1.5	4/0	75 cm	<i>violet</i>	13151524
	DS 24 mm 	2	3/0	75 cm	<i>violet</i>	13201524
		2	3/0	90 cm	<i>violet</i>	13201524-90
3		2/0	75 cm	<i>violet</i>	13301524	
3		2/0	75 cm	<i>violet</i>	13301524-90	
3.5		0	75 cm	<i>violet</i>	13351524	
DS 26 mm 	2	3/0	75 cm	<i>violet</i>	13201526	
	3	2/0	75 cm	<i>violet</i>	13301526	
DS 30 mm 	3	2/0	75 cm	<i>violet</i>	13301530	
	3.5	0	75 cm	<i>violet</i>	13351530	
DS 40 mm 	4	1	90 cm	<i>violet</i>	13401540	
	5	2	90 cm	<i>violet</i>	13501540	

 Reverse cutting		Metric	USP	Length		Order number
	HS 36 mm 	3	2/0	75 cm	<i>violet</i>	13301436
		3.5	0	75 cm	<i>violet</i>	13351436
		4	1	75 cm	<i>violet</i>	13401436

 Blunt loop		Metric	USP	Length		Order number
	HRB 50 mm 	4	1	150 cm	<i>violet</i>	13401850

IN CASSETTES

Use within 4 months after first opening

Metric	USP	Length	Order number
1.5	4/0	15 m	13151
2	3/0	15 m	13201
2	3/0	25 m	13202
3	2/0	15 m	13301
3	2/0	25 m	13302
3.5	0	15 m	13351
3.5	0	25 m	13352
4	1	15 m	13401
4	1	25 m	13402
4	1	50 m	13405
5	2	15 m	13501
5	2	25 m	13502
5	2	50 m	13505
5	2	100 m	135010
6	3&4	15 m	13601
6	3&4	25 m	13602
6	3&4	50 m	13605



PRAXISDIENST
Medical Supplies since 1953

Order here!