

InSight V-IA[®] PLUS

Veterinary Immunoassay Analyser

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes.



PRAXISDIENST
Medical Supplies since 1953
[Order here!](#)



LOW COST,
REFERENCE LAB QUALITY,
IMMUNOASSAY TESTING.

"PASSIONATE ABOUT LABORATORY DIAGNOSTICS"

WOODLEY
EQUIPMENT COMPANY LTD.

InSight V-IA^{PLUS} Veterinary Immunoassay Analyser

An easy to use Immunoassay Analyser providing accurate and reliable results in 3-15 minutes

FEATURES

- Immunofluorescence technology
- Testing time 3-15 minutes (depending on analyte)
- Easy to use
- Accurate and reliable
- Multi-language operation system
- Wi-Fi connectivity
- Reagents stored at room temperature
- Touch screen
- Maintenance free
- Results transmitted via InSight p-LINK to Practice Management Systems



Software Languages

English, Chinese, Spanish, French, Russian, Portuguese, Turkish, Italian, Polish, German

SPECIFICATIONS

Dimensions: 213 (W) x 243 (D) x 195 (H) mm
Weight: 4kg
Power Adapter: 100~ 240V, 50~ 60Hz

Canine Tests Available:

- cPL
- T4
- TSH
- Cortisol
- Progesterone
- NT-proBNP/Troponin I (combined test)
- Cystatin C (SDMA alternative)
- **Core Antibody:** Distemper, Parvovirus, Distemper/ Parvovirus/Infectious Hepatitis Virus combined test
- **Infectious Disease:** Coronavirus, Distemper, Parvovirus, Parvovirus/ Coronavirus combined test
- cCRP
- D-Dimer
- HbA1C
- NT-proBNP
- Troponin I

Feline Tests Available:

- fPL
- T4
- TSH
- Cortisol
- Cystatin C (SDMA alternative)
- **Core Antibody:** Calicivirus, Herpes, Panleukopenia Virus, Herpes/Calicivirus Panleukopenia Virus combined test
- **Infectious Disease:** Coronavirus, Calicivirus, Herpes, Panleukopenia Virus
- D-Dimer
- fSAA
- NT-proBNP
- Troponin I

Equine Tests Available:

- eSAA

Tests In Development:

- fTT4 (Feline Total Thyroxine) Quantitative Test
- HP (Helicobacter Pylori) Quantitative Test
- CPSE
- TOXD Ab
- RB virus Ab
- AFP (Alpha-fetoprotein)
- Fructosamine
- Equine ACTH
- Equine Progesterone
- Foal IgG



+44 (0) 1204 669033 (Option 1)



sales@woodleyequipment.com



www.woodleyequipment.com



Old Station Park Buildings
St. John Street
Horwich
Bolton
Lancashire
BL6 7NY
UK

**"PASSIONATE ABOUT
LABORATORY DIAGNOSTICS"**

**WOODLEY
EQUIPMENT COMPANY LTD.**