

LATEX EXAMINATION GLOVES POWDER-FREE



DESCRIPTION

Teqler latex **examination gloves powder-free** have a wide range of applications for any medical facility, such as medical practices and hospitals. The high quality gloves provide utmost wearer comfort and are skin-friendly thanks to their particularly low protein content (10µg/g).

PRODUCT DETAILS

- High quality latex examination gloves
- Powder-free
- Non-sterile
- Left or right hand, rolled edge
- Tear-resistant
- Micro-beaded surface
- Excellent fit
- High wearer comfort
- Low 10 µg/g protein content, low risk of allergic reaction
- Complies with the European Norm 455, AQL 1.5
- 1 package contains 100 gloves

Overview

Article Number:
T321201 - T321204

Category:
Gloves



Teqler is a trademark of NetMed S.à.r.l.

📍 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

☎ +352 267149 09 📠 +352 267149 19 ✉ info@netmed.lu

Social Media

