


SURGICAL MASK FOR TATTOO ARTISTS



DESCRIPTION

This **surgical mask from Teqler** is particularly suited for tattoo and piercing studios, thanks to the black colour and the dispenser packaging. The surgical mask is optionally available with earloops or ties.

PRODUCT DETAILS

- Black surgical mask
- Compliant with Type II EN 14683
- Bacteria filtration capacity: 98%
- Made of 3-ply non-woven material
- Fibreglass-free
- Soft quality provides high comfort
- With long, flexible nosepiece
- Available with elastic earloops or with ties
- Dispenser box with tattoo design
- Delivery contents: 1 box with 50 masks

Overview

Article Number:
T490112-T329102

Category:
Surgical Masks



Teqler is a trademark of NetMed S.à.r.l.

📍 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

☎ +352 267149 09 📠 +352 267149 19 ✉ info@netmed.lu

Social Media

