

Ambu® Disposable Open Cuff Face Mask



Key Benefits

- Transparent for visual check of bleeding, vomitus, and spontaneous breathing
- Available in 7 different colorcoded sizes for easy size identification
- Soft cuffs and shaped for a comfortable fit for both adults and children
- The masks are phthalate, pvc free, and are not made with natural rubber latex

Ambu® Disposable Open Cuff Face Mask

Ambu® Disposable Open Cuff Face Masks are designed for use with manual and automatic resuscitators and ventilators. Transparent domes for visual check of bleeding, vomitus, and spontaneous breathing. The masks are made with soft cuffs and shaped for a comfortable fit for both adults and children. The masks are phthalate, pvc free, and are not made with natural rubber latex.

Ordering specifications

Item No.	Type
306 1000 00	Ambu® Disposable Open Cuff Face Mask size 1 Neonate
306 2000 00	Ambu® Disposable Open Cuff Face Mask size 2 Infant
306 3000 00	Ambu® Disposable Open Cuff Face Mask size 3 Pediatric
306 4000 00	Ambu® Disposable Open Cuff Face Mask size 4 Small Adult
306 5000 00	Ambu® Disposable Open Cuff Face Mask size 5 Medium Adult
306 6000 00	Ambu® Disposable Open Cuff Face Mask size 6 Large Adult
306 7000 00	Ambu® Disposable Open Cuff Face Mask size 7 Extra Large Adult



Ambu® Disposable Open Cuff Face Mask

Materials

Product	Material
Cuff	TPE
Dome	PolyPropylene



Ambu® Disposable Open Cuff Face Mask sizes 3, 4, 5, 6

Size

Type	Weight
1 Neonate	7.1 g
2 Infant	8.9 g
3 Pediatric	15.1 g
4 Small Adult	18.5 g
5 Medium Adult	21.6 g
6 Large Adult	25 g
7 Extra Large Adult	27.8 g



Ambu® Disposable Open Cuff Face Masks sizes 1-7

Ambu Ltd.
 8 Burrell Road
 St Ives
 Cambridgeshire
 PE27 3LE
 Tel 01480 498403
 Fax 01480 498405
 uksales@ambu.com
 www.ambu.co.uk

Manufactured by:
Ningbo Shengyurul Medical Appliances Co., Ltd.
 No.138 Binhaisi Rd.,
 Hangzhou Bay New Zone
 315336 Ningbo City,
 Zhejiang Province, China

CE
 0120

US: Rx only