

# B. Protected B. Braun

Medibox® Sharps Disposal Containers



CoE Practitioner Use





### Day-to-day safety without compromises

Medibox® sharps containers are designed to efficiently protect health care workers from accidental injuries during sharps disposal.

Together with the B. Braun safety devices, the Medibox® containers help minimizing the risk of transmission of pathogens and contribute to a safe working environment.

- Easy and safe medical sharps disposal
- No-touch, twist-off and insertion inlets for all kinds of used needles, cannulae, scalpel blades and medical sharps
- Handling instructions by self-explanatory pictograms
- Overfill warning by maximum fill line (at approx. 80% of total volume) and visible content level
- Impact and puncture resistant
- Freestanding, ergonomic design
- Easy to operate final lock mechanism



product	sales unit	REF
<b>Medibox®</b>		
<b>Sharps Disposal Containers</b>		
Complying with international standards		
<b>0.7 L</b>	10	9193405
<b>3 L</b>	10	9193413
<b>5 L</b>	10	9193421

## B. Braun Infection Control



Risks of pathogen transmission can never be totally eliminated. Hygienic safety means reducing the number of pathogens, toxic residues and particles to the lowest possible level in order to minimize the risk of nosocomial infections and ensure occupational safety and health.

Due to antibiotic resistances in bacteria and rapid genetic changes in viruses, hygiene products and protective equipment play a crucial role in successfully combatting resistant and pandemic microorganisms.

B. Braun infection control products and services are effectively contributing to the prevention and management of infections in healthcare settings all over the world. Protective wear, hand and skin hygiene, cleaning and disinfection of surfaces and instruments are helping to protect health care workers and patients against all kinds of infectious diseases and to minimize spreading of pathogens.