

ENT

OTOSCREEN 2

Wireless video otoscope



Eclairage Lighting
4 LEDs

1 YEAR WARRANTY



CE



Colour LCD screen



Low battery indicator



Camera



On/off + image capture function switch



Video outlet socket



Protective transport cover

Reference	Model
570 701	Otoscreen 2 video otoscope

References	Descriptions
570 798	Single-use rigid ear speculums, 5 mm (box of 250)
570 794	7.5 cm long camera
570 795	15 cm long flexible camera
570 796	30 cm long flexible camera
570 797	Computer Kit

Colour LCD screen
Integrated high resolution screen (320 x 240 pixels) measuring 60 mm (diagonal) for easier, more comfortable examinations.

Low battery indicator
The indicator warns you when the batteries need to be changed. Otoscreen 2 can run for up to 4 hours with new batteries.

Camera
Otoscreen 2 has a miniature camera allowing you to get closer to the examination zone. Strong lighting is provided by 4 LED lamps. The camera is removable. This means you can take advantage of future improvements to your Otoscreen 2.

On/off + image capture function switch
The same button is used to switch the Otoscreen 2 on and off or to freeze the image. A short press on the button freezes the image. Press it a second time and the device returns to examination mode.

Video outlet socket
Connect Otoscreen 2 on your television with the RCA video wire (included with your item) and check what you see on your screen.

Protective transport cover
Very useful to transport and protect your Otoscreen 2.

Battery compartment
Otoscreen 2 operates with 2 standard LR06 AA 1.5 V batteries which are very easy to change.



Spengler S.A.S. Z.I. La Limoise Rue Robert Maréchal 36100 ISSOUDUN - FRANCE	Tel. + 33 (0) 4 42 52 89 80 Fax + 33 (0) 4 42 52 89 09 export@spengler.fr www.spengler.fr
--	--